



# Edhya TOP LYTE

Equine Dietary feed supplement - Paste form

## ELECTROLYTES - RECOVERY SUPPORT

Edhya TOP LYTE helps compensate the loss of electrolytes in case of important perspiration.

Edhya TOP LYTE also contains methionine, known for its cholagogic action, favouring the use of lipids and dissolution of liver fats.

Vitamin C and E and spirulina input helps muscular recovery and supports muscular tonicity.

### PURPOSE OF USE

- Loss of electrolytes
- Effort recovery
- Transport
- Intensive effort
- High temperatures

**FEEDING INSTRUCTIONS** - Administer directly in mouth.

- 1 syringe (50 ml) after an intensive effort, for 1 to 3 days
- Competition day: 1 syringe on the evening before, and 1 syringe after

Make sure water is available at all times

**Imperative:** ensure the horse drinks water after each effort to control salt concentration in the organism and allow assimilation (a supply on a non-hydrated horse would be pointless)

Before use or before extending the period of use, it is recommended to seek the advice of your veterinary practice.



ELIMINATION  
& RECOVERY



V181220

### COMPOSITION

Dextrose, Sodium chloride, Cider vinegar, Calcium propionate, Lactose, Sodium citrate, Potassium chloride, Sorbitol, Magnesium sulphate, Potassium phosphate, Spirulina

### ANALYTICAL CONSTITUENTS

Moisture	42 %
Crude Protein	3,9 %
Fibre	<1 %
Ash	16,5 %
Chlorine	6,9 %
Sodium	4,4 %
Potassium	1,5 %
Magnesium	0,1 %
Calcium	0,5 %
Glucose	15 %

### ADDITIVES (Per Litre)

#### Amino Acids

DL-Methionine (3c301)	16 g
Glycine	16 g

### Vitamins

Vitamin C (3a300)	5 g
Vitamin E (3a700)	2 000 mg

### Preservatives

Potassium sorbate (E202)	5 g
--------------------------	-----

### Thickening agents

Xanthan gum (E415)	5 g
--------------------	-----

### STORAGE

Store in a cool dry place away from direct sunlight.

### PACKAGING

60 ml syringe  
Box of 10 syringes (60 ml each)



*This is not a veterinary medicine which is subject to authorisation by HPR & IMB*